

Datasheet for ABIN1882417

**AnaTag™ Biotin Protein Labeling Kit**[Go to Product page](#)**1** Publication

## Overview

|              |                |
|--------------|----------------|
| Quantity:    | 1 kit          |
| Target:      | Biotin         |
| Conjugate:   | Biotin         |
| Application: | Labeling (Lbl) |

## Product Details

|                   |              |
|-------------------|--------------|
| Detection Method: | Fluorometric |
|-------------------|--------------|

## Target Details

|              |                                                                                                                                       |
|--------------|---------------------------------------------------------------------------------------------------------------------------------------|
| Target:      | Biotin                                                                                                                                |
| Abstract:    | <a href="#">Biotin Products</a>                                                                                                       |
| Target Type: | Chemical                                                                                                                              |
| Background:  | Biotin-X, SE is amino-reactive biotinylating reagent for peptides and proteins , it has a better avidin-binding affinity than biotin. |

## Application Details

|               |                       |
|---------------|-----------------------|
| Comment:      | Labeling kit          |
| Restrictions: | For Research Use only |

## Handling

|          |      |
|----------|------|
| Storage: | 4 °C |
|----------|------|

## Publications

---

Product cited in: Merchant, Korbelik: "Heat shock protein 70 is acute phase reactant: response elicited by tumor treatment with photodynamic therapy." in: **Cell stress & chaperones**, Vol. 16, Issue 2, pp. 153-62, (2011) ([PubMed](#)).